

Cleaning Wipes Xtreme

Powerful wipes to gently scrub and clean for daily use in the workplace.

We declare war on dirt! Our high-performance fischer Cleaning Wipes Xtreme master every challenge perfectly and remove even stubborn dirt thanks to their practical 2-in-1 function: rough side to loosen the stain and smooth side to remove the dirt.

Product description

High performance hand-cleaning system with a rough and smooth side for quick and skin-friendly cleaning of uncured paints, resins, adhesive, lubricants, oils, tar, inks, coal, graphite, and many other stubborn stains. The cleaning wipes are enriched with propylene glycol and glycerine which work as moisturizers. Dermatologically tested for skin compatibility.


Directions for use

- Before use, remove the lid and peel off the protective foil.
- Pull out a piece of the innermost wipe from the roll and thread it through the hole in the lid.
- Replace the lid and open it.
- Pull up the cleaning wipe (tear it off at the perforation).
- Clean (always test on a hidden spot before use).
- Reseal the tub tightly to maintain moisture and keep it stored upright at room temperature.
- Dispose in the garbage can.

For best results, clean spills immediately. For stubborn stains, soak for five seconds with the wipe. Always test on a hidden spot before use.

Product advantages

- Strong, powerful and grease-dissolving cleaning wipes.
- Reliable cleaning of stubborn dirt.
- Practical 2-in-1 function: rough side to loosen the stain and smooth side to remove the dirt.
- Time-saving thanks to quick cleaning without water and soap.
- Safe and gentle on hands, tough on dirt.
- Dermatologically tested.
- Easy to use thanks to the practical dispenser box with single wipe dispenser.



Properties	Value	Unit	Method / Remarks
	System		
Item Code	567456		Languages on label: EN, DE, FR, ES, IT, PT, NL, NO, SE, DK, FI
Color wipe	Orange		
Odor	Orange		
Grammage wipe	45	g/m ²	
Dimensions wipe	200x200	mm	
Dimensions tub (lwh)	110x110x245	mm	
Dimensions carton box (lwh)	333x223x242	mm	
Content/ tub (wipes)	100	pcs	
Content/ carton box (tubs)	6	pcs	
Weight tub	714	g	
Weight carton box (net)	4,284	g	
Shelf life unopened	31	months	
Shelf life opened	6	months	
Dermatological Test	Fresenius Institute		
pH value	5.5-6.5		

Warning

If medical advice is needed, have product container or label at hand. Keep out of reach of children. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: Get medical advice/attention.

For further information on handling, please refer to the safety data sheet; this contains important information.

The application instructions in this technical data sheet and our technical advice on application, whether spoken or written, are given to the best of our knowledge but do not constitute either a warranty as to the nature or usability of the products or an independent or dependent promise or warranty declaration of any kind. We recommend to always check the suitability of our products for the intended use. Due to the wide variety of possible applications for each individual product and the unpredictable conditions at the place of processing, we recommend testing the function on the substrates before use.